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Conductimétrie

Chemical measurement

Conductimetry

Réf :
701 330

Français – p 1

English – p 5

Version : 7006

Sonde conductimétrique CTA

CTA conductimeter probe

1 Description

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0 825 563 563*

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1 Description

The CTA conductimeter probe equipped with graphite electrodes is robust and easy to maintain.

It is equipped with a temperature sensor for allowing a usage in automatic temperature compensation with the conductimeter CTA Initio® Ref. 701 328.

Diagram 1

2 Technical characteristic

Electrodes	Graphite
Cell constant	Theoretical ⁽³⁾ : 0,393 cm Experimental ⁽⁴⁾ : Measured at 5kHz : 0,64 - 0,69 cm
Minimum height of liquid required	20 mm
Dimensions (mm)	∅ body: 12; ∅ tip: 16 Length: 140
Temperature sensor	CTN Probe
Connection	Mini DIN socket
Measurement range	0 - 20 mS.cm ⁻¹ (Connected to CTA Initio® conductimeter)

2.1 Geometry of the probe


Diagram 2

⁽³⁾ Calculated on the basis of relation $K = S/L$. S is the surface of an electrode and L is the distance between the electrodes.

⁽⁴⁾ Average values measured on a batch of 10 probes at 5 kHz, frequency of operation of the CT. conductimeter

2.2 Resistance to aggressive chemical products

This probe is designed for resisting the immersion in saline, acid, basic, diluted electrolytic solutions.

Do not dip this probe in oxidising solutions (potassium permanganate, sodium hypochloride, etc.) or in organic solvents.

For an optimum service life, it is recommended to always wash the probe with distilled water after each usage.

3 Connexion

This probe is equipped with a cord provided with mini DIN socket.
It is connected to the conductimeter CTA (Ref.: 701 328).

4 Stability of measurements

If the probe is not used for a long time and/or if it is stored in a dry ambience, it is recommended, before a series of measurements, to humidify the surfaces of graphite by dipping the electrodes in distilled water for a few minutes (5 to 30 minutes depending on the degree of drying). This ensures the stability of measurements obtained.

5 Maintenance

Rinse the tip of the probe with demineralised water after each usage (do not use the organic solvents).

6 After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.
For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
468 Rue Jacques Monod
CS 21900
27019 EVREUX CEDEX France

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediately to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.

** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux