Physiologie

Physiology


**Ref :
512 021
512 027**

Français – p 1

English – p 9

Version : 6010

**Maquette fumeur passif
*Passive smoker model***

1 Introduction

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
NOTES

1 Introduction

The purpose of the model is to make students aware of the health hazards posed by tobacco use to the smoker himself as well as to the people around him. The “passive smoker” model only works simultaneously with a “mechanical active smoker model.”

The harmful effects of tobacco on the lungs are illustrated on a piece of filter paper placed in the “chest cavity.”

As the smoke is being “inhaled” or aspirated, one can directly observe the tar being deposited as the smoke passes through this paper.

2 How it works

Have the “active smoker” smoke. The air contained in the enclosure will be filled with smoke coming from the burning of the cigarette as well as the smoker’s exhalations.

Now have the passive smoker breathe.

In both cases, the air or the smoke from the cigarette is “inhaled” with the aid of a syringe. Before entering the syringe, the smoke-laden air passes through a double compartment. This double chamber is divided by a partition consisting of a sheet of filter paper inserted between two discs. The discs have a central cutout, the outline of which is in the shape of the lungs. The smoke inhaled must pass through the paper via this cutout. The tar is retained on the paper and therefore causes it to turn brown. The zone thus stained will have the silhouette of the lungs.

Note: This device functions equally well with tobacco-free and nicotine-free cigarettes.

3 Contents

The **kit consisting of the active and passive smoker models (order no. 512 027)** contains the following:

- an **active smoker model** (see the description below)
- a **passive smoker model** (see the description below)

3.1 Active smoker

The **active smoker model (which can also be ordered separately under item no. 512 045)** consists of the following elements (see Figure 1):

- a 46 x 20 cm support panel silkscreened on both sides **(5)**,
- two 30 x 6 cm rectangular plates with cutouts that represent the feet **(5)**,
- a transparent chest cavity for the assembly of:
 - The front cover of the lung chamber, equipped with a single orifice **(2)**,
 - Two lung discs 70 mm in diameter **(3)**,
 - The back cover of the lung chamber, equipped with two orifices with different diameters **(4)**,
- a 60 ml plastic syringe **(6)**,
- 3 connecting tubes:
 - one flexible tube with a large diameter,
 - two flexible tubes with small diameters.
- 10 filter paper discs 70 mm in diameter **(7)**.

3.2 Passive smoker

The **passive smoker model (which can also be ordered separately under item no. 512 021)** consists of the following elements (see Figure 2):

- a transparent enclosure that is foldable with a silkscreening of a “passive smoker” **(5’)**,
- a transparent chest cavity for the assembly of:
 - The front cover of the lung chamber equipped with a single orifice **(2)**,
 - Two lung discs 70 mm in diameter **(3)**,
 - The back cover of the lung chamber, equipped with two orifices with different diameters **(4)**,
- a 60 ml plastic syringe **(6)**,
- connecting tubes:
 - one flexible tube with a large diameter,
 - two flexible tubes with small diameters.
- 10 filter paper discs 70 mm in diameter **(7)**.

4 Assembly

4.1 Assembling the active smoker

Attention: The assembly and the connections of the active smoker model are slightly different when using it together with the passive smoker model and when using it by itself without the passive smoker enclosure (for the latter use, consult instruction manual 512 045).

- First, locate the front **(2)** and back **(4)** covers of the lung chambers (see Figure 1).
 - Attach the support panel to the two feet, in the slots, as shown in Figure 1.
 - On the front side of the support panel (the side with the human model), press the front cover of the lung chamber **(2)** into the round opening provided. The alignment slot and the orifice are oriented towards the top (see Figure 1).
 - Then make the following assembly in the back cover **(4)**:
 - Insert a lung disc **(3)** in the cavity of the back cover **(4)**. An alignment slot only allows one position for the disc.
 - Place a piece of filter paper **(7)** on the lung disc and then cover it with the second lung disc **(3)**.
 - Slide this assembly into the front cover of the chamber attached in the preceding step; it should protrude from the back side of the panel (the side silkscreened with "the dangers of tobacco"). The lungs are now in place.
 - Connect the flexible tubes, feeding them through the openings provided on the support panel (see Figure 1).
 - Connect the syringe **(6)** to the flexible tube coming from the **left hip of the model**. The syringe is operated from **the back of the active smoker** (outside the passive smoker enclosure).
 - Feed the flexible tube **(A)** through the slot provided by the left knee; it is operated from outside of the enclosure (see Figure 7) and is designed to recycle the contaminated air to the front of the model and thus into the interior of the enclosure.
 - Place a cigarette **(1)** in the end of the flexible tube representing the mouth.
- The model is now ready to be outfitted with the passive smoker enclosure and to operate.


Figure 1

4.2 Assembling the passive smoker

- Assemble the enclosure as shown in Figure 2: fold the four sides of the enclosure and insert the flaps in the slots provided. Leave the lid open in order to continue with the assembly.
- Note: It may be easier to perform the following assemblies with the enclosure unfolded: In this case, locate the future inside and outside faces of the enclosure, assemble the lung chamber, connect the tubes and complete the assembly by closing the enclosure again.
- Press (from the inside face) the front cover of the lung chamber **(2)** in the round opening on the support panel. The alignment slot and the orifice are oriented to the top (see Figure 3).
- Then make the following assembly in the back cover **(4)**:
- Insert a lung disc **(3)** in the cavity of the back cover **(4)**. An alignment slot only allows one position for the disc.
- Place a piece of filter paper **(7)** on the lung disc and then cover it with the second lung disc **(3)**.
- Slide this assembly into the front cover of the enclosure; it should protrude from the outside face of the enclosure. The lungs are now in place (see Figure 2).
- Connect the flexible tubes, feeding them through the openings provided on the enclosure (see Figures 4 and 5).
- Connect the syringe to the flexible tube exiting the enclosure near the left hip of the model (see Figure 6).
- The syringe must be placed outside the enclosure so that it can be operated.
- Feed the flexible tube **(A')** through the slot provided by the left knee; it is operated from the outside of the enclosure (see Figure 6) and is designed to recycle the contaminated air to the front of the model and thus into the interior of the enclosure.
- Close the enclosure again.
- The model is now ready to be attached to the active smoker model (see Figures 6 and 7).

Assembly (outside view)


Figure 2

Assembly (inside view)


Figure 3


Figure 4

Complete assembly; view of the passive smoker


Figure 5

Complete assembly; transparent view of the active smoker


Figure 6


Figure 7

5 Operation

Set up the smoker and passive smoker assembly near a window that can be opened all of the way.

Place a protective sheet of paper on a tray or on the work surface to catch the ashes falling from the burning cigarette.

You can make the active smoker smoke and the passive smoker breathe simultaneously.

5.1 To make the active smoker smoke

- Place a cigarette in the end of the flexible tube representing the mouth and light the cigarette.
- Close the cover of the enclosure. Because the enclosure is never completely sealed (it does limit the amount of smoke emitted into the classroom), it is possible that some smoke may escape from the top; however, most of the smoke should remain confined.
- The plunger of the syringe should be in the depressed position.
- Pinch the flexible tube **(A)**.
- Pull on the plunger of the syringe. The aspirated smoke will pass through the filter, and you can observe it turning brown.
- Release the pressure on the flexible tube **(A)** and depress the plunger again. The contaminated air will now be expelled via the tube A and recycled into the enclosure.
- Repeat the two preceding operations.
- These operations can be repeated until the cigarette is consumed.
- The air in the enclosure becomes filled with smoke.

5.2 Making the passive smoker breathe

- The plunger of the syringe should be in the depressed position.
- Pinch the flexible tube **(A')** entering the enclosure at the left knee of the model (see Figure 6).
- Pull on the plunger of the syringe. The aspirated smoke will pass through the filter. You can observe it turning brown.
- Release the pressure on the flexible tube **(A')** and depress the plunger again. The contaminated air is recycled into the enclosure.
- Repeat the two preceding operations.

Once the demonstration is completed and the cigarette is consumed, the smoke-filled air should be emptied from the enclosure.

Open the window to avoid filling the room with smoke; quickly lift up the active and passive smoker assembly (by holding it by the feet and the enclosure) and quickly point the open base towards the window to let the contaminated air escape outside the classroom.

6 Teaching and technical suggestions:

- The face to face positions of the active smoker and the passive smoker can simulate two individuals having a conversation: the one smoking and the other not smoking and thus believing that he is safe from the damage caused by tobacco. This example can form the basis for a classroom discussion.
- Also, the confined enclosure and the view of the smoke accumulating therein can lead to a discussion of tobacco use indoors or in public places, etc.
- If you have a high precision scale, you can measure the following:
 - Weigh the filter prior to assembly, repeat the aspiration until the cigarette is entirely consumed.
 - Weigh the contaminated filter. Its increase in mass is a measure of the tar deposited. Compare this with the value given on the pack of cigarettes.

Note: This measurement is skewed by the weight of the combustion water retained by condensation on the filter.

7 Storage and care

After use, take the lung chambers apart and clean them with denatured methanol, and then with water. Discard the contaminated filters. Clean the enclosure with a moist sponge. If necessary, use the scouring side with care to remove any smoke residue.

8 Consummable materials

Pack of 100 filters with a diameter of 70 mm

Item no. 703 296.

9 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediately to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service * :
+ 33 (0)2 32 29 40 50

** Call cost not included*