Microscopes

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Model shown : LM900 Wide field Ref : 571142

1 Description

	Concerned models	
	PCB Type Ref : 571 141	LM 900 Ref : 571 142
Robust and stable stand inclinable through 90°	☒	☒
Round stage plate Ø 120 mm, rotatable through 360°, centrable with fixing clamp	☒	☒
Stage centering device consisting of 2 side knobs	☒	☒
Revolving turret for 3 objectives	☒	☒
Optical tube length of 170 mm	☒	☒
Huyghens eyepieces X 6 and X 15 and wide fields x 10	☒	☒
Achromatic objectives with international diameter of 20 mm 4 x (NA 0.10) 10 x (NA 0.25) 40 x retractable (NA 0.65)	☒	
Achromatic objectives with international diameter of 20 mm 4 x (NA 0.10) 10 x (NA 0.25) 60 x retractable (NA 0.85)		☒
Focussing through rack and micrometer screw	☒	☒
Reflector : plane and concave mirror	☒	☒
Iris diaphragm	☒	☒
2 lens Abbe condenser (NA 1.2)		☒
Filter holder + blue disk		☒
Weight of microscope	4.000 kg	4.300 Kg
The microscopes are delivered in polystyrene boxes with the manual		

2 Technical features

- Magnifications

Huyghens eyepiece "H" Wide field « WF »	Achromatic objectives		
	4 X	10 X	40 X ou 60 X
6 x H	24	60	240 or 360
10 x WF	40	100	400 or 600
15 x H	60	150	600 or 900

3 Operating instructions

3.1 Lighting device

The PCB microscopes are fitted with a plane/concave mirror with harp, adjustable on two perpendicular planes.

The light source could be :

- Either natural (day light but preferably not direct rays from the sun).
- or artificial (for example an individual lamp : Ref: 574 169).

For longer viewing, we recommend the use of the plug-in lamp instead of the mirror (description given in "ACCESSORIES").

3.2 Observation in natural light (mirror or lamp)

The sample, placed on the stage plate, is held by two clamps.
Start the observation with an objective of small magnification (4 x or 10 x)
Direct the light beam on the sample by adjusting the direction of the mirror.
Vary the light intensity by adjusting the diaphragm so as to obtain the best contrast.

Do a rapid focus by using the adjusting screw of the rack. If required, bring the focus to perfection with the help of the micrometer screw. To move to a higher magnification, just turn the turret to bring the desired objective above the sample.

Focus by only using the micrometer screw .

Note : The objectives 4 x and 10 x can be used for observing a sample without cover glass.

These objectives are also suitable for observation in reflected light.

Lighting conditions depend greatly on the position of the condenser : remember to adjust it without the sample to obtain a uniform field.

4 Accessories (not delivered with the microscopes)

4.1 Lighting device

Plug-in lamp with earth plug (Ref. 574 169) 220 Volts

4.2 Observation in polarised light

For more accurate observations use the POLARISING STAGE, (Ref. 574 174) which will give you best results.

5 Observation with cameras

Use the wide field eyepiece x 10 delivered with the microscope.
For more details refer manual delivered with the cameras.

6 Spare parts

Mirror with harp (Ref. 574 190)
Eyepieces and objectives: Refer JEULIN catalogue

7 Maintenance

The microscope must be kept away from dust and from possible shocks.
When not in use, keep it in the box or cover it with the help of a dust cover.

The metal parts of the microscope must be cleaned regularly with a soft cloth soaked in acid free Vaseline. If required, clean with spirit before lubricating.
During cleaning and lubricating, take care not to dirty the optical parts of the microscope.

In order to ensure a dust free prism, never remove the eyepieces.

Never touch the lenses.

If dust has entered the objective, return it to us for cleaning.

Before transport, pack the microscope and its accessories.

8 After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.

For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediately to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service * :
+ 33 (0)2 32 29 40 50

** Call cost not included*