TP de chimie

Chemistry tutorial

Ref :
702 082

Français – p 1

English – p 3

Version : 6010

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This set of three kits is designed for the teacher. It can be used for representing the close-packed or split model of different molecules and for solving the problems in stereochemistry.

1 Composition

Qty.	Element		Colour	No. of holes / type	Ø mm
30 12	Carbon	C	Black	4 / tetrahedron sp ₃ 5 / dipyramid	23
60 4	Hydrogen	H	White	1 2 / linear	
6 6 4	Oxygen	O	Red	4 / tetrahedron 2 / angular 1/	23
8 4	Nitrogen	N	Blue	4 / tetrahedron 5 / tripyramid	23
4 8	Sulphur	S	Yellow	4 / tetrahedron 2 / angular	23
8	Chlorine		Green	1	
2	Bromine		Orange	1	
1	Iodine		Purple	1	
3 2	Phosphorous	P	Mauve	4 / tetrahedron 5 / tripyramid	23
3	Metal		Grey	6	
1	Alkyne group pre-mounted				
1	Alkyne group pre-mounted				
40	Medium link		Grey		19/31
12	Long flexible link		Grey		32/43
60	Link (close-packed model)		Transparent		
20	"Accordéon" link		Grey		
60	Short link		White		2/10
1	Dismantling tool				

2 Types of links

- The grey links are used for symbolising the simple covalent links in split model. The scale of representation of inter-nuclear distances in the split model is approximately 2.5 cm/Ångström.
- The long flexible grey links allow representing the multiple covalent links in the split models (double or triple links).
- The short white links are used for associating the hydrogen atoms 1 hole with other atoms. They can however be used between any other type of atoms instead of the medium size links for developing close-packed models. The scale of representation of inter-nuclear distances in the split model is approximately 1.5 cm/Ångström.
- The transparent links allow associating any type of atoms, except the hydrogen atoms 1 hole, in close-packed model. The scale of representation of inter-nuclear distances in the split model is then approximately 1.3 cm/Ångström.

- The grey "accordéon" links are used only for building cycles in which the atoms are coplanar

3 Using the dismantling tool

- Place the tool on the link to be removed (the face with the markings should be turned upwards).
- Push the key horizontally and completely under the link which will be lifted by 2 mm. Do not try to lift the tool before positioning it correctly; this will create unnecessary tensions on the edges of the dismantling tool.
- Release the link gently through a lifting movement on the tool while holding it with the thumb.

Note:


The scale of representation is not respected in the iodine and bromine atoms that are of small size. They are used for studying the stereochemistry (chirality).

4 After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.

For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50

*The Molymod[®] molecular models are made in England by SPIRING Enterprises Ltd.+
The registered trademarks MOLYMOD[®], MOLYDOME[™] and ATOMOD[™] are the exclusive properties of Spiring Enterprises Ltd and all rights are reserved.*

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediately to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service * :
+ 33 (0)2 32 29 40 50

** Call cost not included*