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Periodic phenomena

**Ref :
222 023**

Français – p 1

English – p 3

Version : 7001

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Vibrator for Melde's experiment

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This is the basic device for the study of transverse and longitudinal standing waves in vibrating motion.

The observation of transverse standing waves along a string constitutes Melde's experiment; the vibrator produces the vibration that is electrically sustained.

We easily obtain longitudinal standing waves with a spring.

All these experiments help verify the phenomenon of high-pitched resonance. Melde's experiment also helps verify the formula for vibrating strings:

$$f = \frac{k}{2l} \sqrt{\frac{F}{\mu}}$$

1 Principle – description

1.1 Principle

The vibrator is essentially made up of a small electromagnet and a strip that vibrates at 100 Hz when excited directly by the mains voltage.

1.2 Accessories

Accessories required for experiments, not supplied with the device:

Mandatory equipment

- 2 rod supports, weighing at least 3 kg (a lighter support will have to be loaded or fixed to the table with a clamp, for example).
- A-shaped stand with rod (Part no.: 701031)
- 2 tightening nuts tightening nut (Part no.: 703099)
- A thin cord of a visible colour 25 metres of thin yellow cord (Part no. 323011)
- 1 Spring for standing waves – diameter 10 mm – length 800 mm (Part no.: 222002)
- 1 graduated metre rule

Useful supplementary equipment

- 1 pulley on a rod (Part no. 323017)
- 1 Box of weights with hooks (Part no. 703016)

Melde's vibrator (Part no. 222023) complies with the EN61010-1 standard installation category II, pollution degree 2.

2 Experiments

2.1 Transverse vibrations of a string

Place two supports at about 1 metre from one another; fix the vibrator on one of them using a nut in such a way that the vibrating strip is horizontal and directed towards the other support. On the latter, fix the pulley using the 2nd nut, at the same level as the vibrator.

Then, attach one end of the string in the hole of the strip, and the other to the load-bearing pan. Pass the string over the pulley.

You will obtain maximum amplitude for a given value of the tension in the string. For example, you obtain 2 anti-nodes with an additional mass of 200 g.

2.2 Longitudinal vibration

Place a support so that its rod is close to the edge of the table. Fix the vibrator at the end of the support rod so that the end of the strip faces the operator and the screen-printed face is upwards.

Hook the spring so that it can hang freely in front of the table.

You could also attach the end of the spring to a fixed point, whose distance determines the tension of the spring.

Adjust to obtain antinodes of maximum amplitude.

3 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service * :
+ 33 (0)2 32 29 40 50

* Call cost not included