Didactique

Kit de chimie

Didactic

Chemistry kit

Ref :
106 074

Français – p 1

English – p 6

Version : 0106

Kit miroir d'argent

Silver Mirror Kit

1 Description

1.1 Généralités

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1 Description

1.1 General

The silver mirror kit allows conducting several experiments in minimum time. With pre-dosed packaging the teacher can devote the entire time to teaching rather than preparation. The kit comprises dropper bottles secured prior to first use by leaktight sealing, thus minimising risks in the laboratory.

The concentration of reagents present in the dropper bottles is reduced to a minimum to eliminate risks when handling the products. All the products in the silver mirror kit are free of organic solvents and do not emit any volatile organic compounds (VOC) so can be used in the classroom and do not require the use of an extractor hood.

1.2 Description of Composition

The silver mirror kit comprises 2 bottles, one of ammonia and the other of silver nitrate and 5 micro-tubes containing glucose.

This kit allows performing the following 5 full experiments:


Photo 1

- 1- Glucose micro-tubes
- 2- 1 glass 30 ml dropper bottle of silver nitrate
- 3- 1 PE 20 ml dropper bottle of ammonia

Table 1

2 Additional Equipment

The equipment necessary for satisfactory production of the silver mirror is as follows:

Designation	Quantity
Test tube 16 or 18 diameter	1
A heating element	1
A small graduated flask or pipette	1
Pair of protective gloves in vinyl or latex ⁽¹⁾	1
Pair of safety glasses ⁽¹⁾	1

Table 2

(1) Individual protection equipment recommended although not mandatory

3 Usage and safety precautions

3.1 Warning

The Jeulin Company cannot be held responsible for an accident that occurs after the iron protection and corrosion kit has been used under conditions other than those described in these instructions.

Similarly, the Jeulin Company cannot be held responsible for an accident that occurs because the safety instructions given in these instructions are not respected.

3.2 Labelling and safety datasheet

Flasks containing the reagent are labelled individually in accordance with the regulations in force.

Before performing any manipulation, get pupils to read the label so that they can be familiar with the risks and protection measures to be applied.

If any additional information is required, the safety datasheet for the solution can be downloaded from the www.jeulin.fr site

Product	FDS reference
Silver mirror kit	106 074

Table 3

FDS = Safety Datasheet

3.3 Individual protection

The iron protection and corrosion kit was formulated to minimise risks related to use of the product. Although not compulsory, we recommend that safety goggles shall be worn in addition to the lab coat during handling. It is also recommended that gloves should be worn.

3.4 Collective protection

No collective protection equipment is required during storage (safety cabinet, etc.) or during use (hood, etc.), because reagents do not contain any flammable volatile or toxic solvents.

3.5 Waste


The iron protection and corrosion kit does not contain any components that create any particular danger for the environment. However, these products must be considered as being chemical laboratory waste and treated as such, in the context of a citizen environmental protection policy. Reagent quantities used in this kit are kept low in order to minimise waste volumes and disposal costs.

4 Procedure


4.1 Principle of the silver mirror reaction

The silver mirror kit allows manufacturing a silver mirror by oxidising silver ions in an aldehyde, in this case glucose. The silver ions are precipitated by heating a solution of ammoniac silver nitrate in the presence of aldehyde. The precipitate attaches to the walls of the glass test tube to form a silver mirror.


Obtaining the silver mirror is based on two principles. First of all it is necessary to manufacture the ammoniac silver nitrate in two stages. Then react the aldehyde with the ammoniac silver nitrate.

Addition of ammonia to silver nitrate. A sodium hydroxide precipitate is formed.

When the precipitate is dissolved, addition of ammonia stops. The ammoniac silver nitrate is formed.

During the reaction, the aldehyde reduces the Ag^+ ions to form solid Ag.

At the time of addition of aldehyde, and under the effect of heat, the solution becomes brown then black. The silver precipitate is formed on the walls of the test tube.

4.2 Preparation of the operation

Assemble at each work station the necessary accessories and consumables in the kit.

Photo 2

4.2.1 Preparation of ammoniac silver nitrate

Cut off the top of the PE dropper bottle containing the ammonia (photo 3). Then pour 5 mL of silver nitrate into the test tube (photo 4). Add to this reagent 3 mL of ammonia. For correct dosage, it is preferable to use a graduated flask. The precipitate of silver hydroxide is obtained (photo 5).

Photo 3

Photo 4

Photo 5

4.2.2 Obtaining the silver mirror

Open a micro-tube of glucose. Pour all the contents into the test tube and mix (photo 6). Light the electric heater element and place the tube close to it. **Direct the tube towards a window or a wall to control projection of the product.** Agitate gently during the heating cycle to avoid heating only a proportion of the reaction solution and caramelising the glucose (photo 7).

Photo 6

Photo 7

4.2.3 Result

After a relatively short heating cycle depending on the adjustment of the heater element, a silver mirror is obtained (photo 8).

Place the test tube on a holder and leave to cool (photo 9). Once cooled, empty the reaction liquid remaining in the tube into an appropriate receptacle to avoid spilling liquid during multiple operations.

Photo 8


Photo 9

5 Cleaning of the equipment

To clean the test tube, fill with water and washing up liquid. Leave to soak for some time before rubbing. Recovery of all the silver from the glass walls is not guaranteed unless the silver precipitate is soaked in nitric acid, even of a low concentration.

6 After-Sales Service

For any repairs, adjustments or spare parts, please contact:

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service *

+ 33 (0)2 32 29 40 50**

* Call cost not included

** Only for call from foreign countries