Electricité

Electrostatique

Electricity

Electrostatics

Ref :
273 004

Français – p 1

English – p 5

Version : 4102

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Electrostatics kit

1 Composition

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1 Composition

The kit is made up of :

- A pendulum mounted on a small frame (gerglas support beam),
- A very light bowl recovered by a conductive layer and suspended by a silk thread,
- A hard rubber stick \varnothing 10 mm (200 mm long),
- A glass stick \varnothing 12 mm (200 mm long),
- A piece of kidskin.
- A piece of synthetic fiber

2 Aim of the kit

- Show evidence of electrifying phenomena by friction.
- Show the existence of '2 electricities', positive and negative.
- Introduction of the notion of electric charges.
- Show evidence of electric forces.

3 Use

3.1 Pendulum set up

Strongly push in the rubber stopper mounted on the rod into the hole in the plate, so that the support beam is found, preferably, parallel to the large side of the plate.

3.2 Characteristics of elements

- The pendulum ball is made up of a shot of graphite cellulose.
- The hard rubber stick must be rubber with the kidskin ; using electron capture pulled towards this, the surface becomes seat for negative charges.
- The glass stick must preferably be rubbed with a wool fabric (cloth or tight knit fabric) From loss of electrons from the wool, its surface becomes a carrier of positive electronic charges.

3.3 Experiments

❖ Experiment I : Showing evidence of electrifying by friction

Energetically rub the hard rubber stick with the kidskin and put it near to a metal object : you can hear the characteristic noise of a small spark which can even be seen in darkness. This spark proves that the rubbing has provoked the appearance of an **electrical** phenomenon in the stick.

Bring the rubbed stick near to pieces of paper. They are attracted: creation of an electric **field** giving rise to electrical **forces**.

This field can be made "physically" seen if you put the electrical stick close to a face: the standing up on end of the facial hair creates a sensation characteristic of a 'spider's web'.

The same experiment can be done with the glass stick which must be **dry**.

❖ Experiment II : repulsion of charges of the same sign

1). Rub the hard rubber stick with the kidskin and gradually approach the pendulum ball : this is attracted (phenomenon already shown in another form in experiment I).

As soon as the ball touches the stick, it is repulsed: on contact with the stick the charges of it are passed to the surface of the ball which is then charged with electricity of the **same sign**.

2). After having discharged the ball by gently touching it between 2 fingers, do the same experiment with the glass stick that has been rubbed with wool cloth : the same phenomena are produced.

❖ Experiment III : distinction of hard rubber and glass electrifying

Charge the pendulum ball as in experiment II (by touching it with the hard rubber stick).

Quickly move the rubber stick away as soon as the ball has been repulsed. Rub the glass stick and move it towards the charged ball : the ball is greatly attracted and does not immediately detach from the glass stick. The charges, carried by the glass exert an **attraction** on the charged ball, are therefore of the opposite sign to those having been captured by the rubber stick.

4 Electronic interpretation of phenomena

The work on friction forces provoke 2 phenomena:

- Heating of the object being rubbed,
- Movement of unstable electrons in the peripheral layers of atoms making up the present materials. These materials being isolated, the excited electrons cannot be spread into the mass of the material and be discharged towards the ground by the electrical chain made up by the operator. They therefore remain confined to the neighbouring surfaces in contact and in their agitation, they leave one surface to go to another: the direction of this movement depends on complex factors such as the density of electron clouds mobilised in each material, their respective molecular structures, the arrangement of atoms in these molecules, etc...

It is always such that a body is 'negatively electrified' which means that it possesses **too many** electrons; it is therefore willing to give them up at the first opportunity to a body with too few electrons.

It is also necessary to point out that a single negative charge has a real physical presence: it is represented by an electron (you can also say a negatron), particle that can be individualised and certain other characteristics we have determined (mass, dimension, etc...).

A positive charge, at the level we are at, does not really exist, it is **only a manifestation of a lack of negative, negatron, charges**, a 'hole to plug'. There are no positively 'charged' bodies but a body that lacks negatrons to find their electrical equilibrium.

5 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediately to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.
** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux