Electricité

Magnétisme

Electricity

Magnetism

Réf :
292 004
292 005
292 006

Français – p 1

English – p 3

Version : 1109

Circuits types

Circuit types

1 Principe – description

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0 825 563 563 *

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
1 Principle – description

The apparatus is designed for observing the magnetic force lines which produce the fields created by a current in an element of the rectangular or circular circuit, or a solenoid.

The direction of the lines of force can be revealed by using small magnetic needles (our reference 263 010).

Assembling circuits on a transparent altuglass base allows the lines of force to be projected onto a screen or white wall. Use an overhead projector or an omniscopes. The back side of the different circuits can be easily removed in order to be used with an overhead projector and to allow a display for all the students.


The complete set consists of three pieces of apparatus:

Projectable solenoid, photo 1

Reference 292 006

18 turns of 10/10 wire enable you, with a 7 A current, to obtain complete magnetic fields which clearly show the internal and external lines of induction from a magnetic field.

Circular current, photo 2

Reference 292 005

The spool consists of 25 turns of relatively fine wire (5/10). It thus forms a core of fairly small section to give students the idea of the "single loop" while avoiding the high currents which would be necessary to obtain a magnetic field with a single wire.

A current of 4 A is enough to get a complete magnetic field.

Rectangular current, photo 3

Reference 292 004

Rectangular frame, of which the small side is used. The winding contains 7 turns of enamelled 10/10 wire, thus with a current of 6 A it gives a circular field. With magnetic needles, a lower current is sufficient.


1.1 Supplementary accessories

- Individual supply
- 6 A rheostat
- Knife switch
- 1 kg of iron filings
- Omniscope
- Overhead projector

2 Experiments




Important

To obtain a good pattern, it is pointless exceeding the stated currents. It is thus necessary to put a rheostat into the supply circuit to control the current. Without this precaution, you can cause the wire to heat which can damage the apparatus.

1. Rheostat
2. Supply
3. Circuit type

Projection requires use of an omniscope or overhead projector.

3 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediately to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service *

+ 33 (0)2 32 29 40 50**

* Call cost not included

** Only for call from foreign countries