Physiologie

Physiology

**Ref :
512 024**

Français – p 1

English – p 3

Version : 6010

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Pulmonary ventilation model

1 Composition

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1 Contents

- 1 transparent, polycarbonate bell jar
- 1 latex membrane and fastener
- 2 rubber balloons
- Joints and couplers

2 Notes

Pulmonary ventilation (breathing) is a purely mechanical process. Air enters the lungs when the diaphragm contracts and causes the rib cage to expand. This action causes the volume of the lungs to expand via the pleura and fill with air. When the diaphragm relaxes, the lungs contract and the air is expelled. During physical exertion, the muscles attached to the ribs also intervene and increase the volume of the rib cage, and consequently the amount of air inhaled, even more.


This model only illustrates the overall mechanism of increasing the volume of the rib cage that results in the expansion of the volume of the lungs. The various diaphragm and rib movement mechanisms are beyond the scope of this model.

3 Instructions

By pulling on the membrane, the volume of the bell jar is increased and the pressure inside the bell jar decreases, which results in an increase of the volume of the balloons and thus their inflation by air being drawn in.

By releasing and by pushing on the membrane, the volume of the bell jar is reduced and the pressure inside the bell jar increases, which results in a decrease of the volume of the balloons and thus their deflation by air being expelled.

These mechanisms allow one to explain what happens in the rib cage and thus the manner in which air is drawn in and expelled from the lungs.

4 Storage and maintenance

- Store in a dark place, as light (direct sunlight, ultraviolet light) may cause the membrane to deteriorate.
- Store in a place protected from heat.

5 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service * :
+ 33 (0)2 32 29 40 50

** Call cost not included*