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Physics PW

Fluids mechanics: JURIN's Law

Réf :
243 015

Français – p 1

English – p 3

Version : 9101

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Capillarity Study Device

1. Objectifs

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
1. Objectives

- Observe the capillarity phenomenon in a perfectly wettable capillary tube.
- Verify Jurin's law and calculate the surface tension coefficient (or the surface tension).

2. Description

This device includes:


- 1 tube holder provided with a tank, entirely in plexiglass
- 4 glass capillary tubes of internal diameter 0,30 mm ; 0,55 mm ; 0,9 mm ; around 1 mm (be careful, it may be subject to manufacturing variations)

A filling tube facilitates the introduction of liquid into the tank. The back of the frosted support enables a backlight to improve reading.

Tubes length:	250 mm
Support dimension:	200 x 116 x 29 mm
Tank dimension:	40 x 116 x 29 mm

3. Principle

When a capillary tube is immersed into a wetting liquid, it rises in the tube and "hangs" to the inner wall of the tube by forming a meniscus.

Jurin's law calculates the ascension h of a ρ density liquid in a capillary tube of internal diameter $2r$:

$$h = \frac{2A}{r \cdot \rho \cdot g} \cdot \cos \alpha$$

A : surface tension in N/m

α : angle of liquid-wall connection

water (perfectly wetting liquid) $\alpha = 0^\circ$

4. Manipulations

4.1 Experimentation

Mount all the capillary tubes on the device support.

Fill the device tank with water.

Note the water temperature.

Observe the water ascension in the capillary tube.

Observe the water-air separating surface in the capillary tube.

To perform measurements, it is necessary to get the capillary tube "wet". For this, draw the liquid back into each capillary tube and let the liquid go down.

4.2 Exploitation

Measure the water ascension height h in each tube.

To improve the reading, it is possible for example to:

- illuminate the back of the device holder;
- have a sheet of graph paper on the support;
- use a colorant.

Prove that the water ascension height h in a capillary tube of internal radius r is proportional to $\frac{1}{r}$.

Calculate the water surface tension A in contact with air at the temperature of the experiment.

Results:

Water temperature:

Internal diameter $2r$ (mm)	0,30	0,55	0,90	1
water ascension height h (mm)				
$h r$ (mm ²)				
A (N/m)				

Perform the same experiment using different liquids.

Values of surface tension for some liquids in contact with air is at 20°C:

Liquid	A (10 ⁻³ N/m)
ether	17
ethanol	22
Teepol solution at 1%	45
glycerin	63
water	72,5

5. After-sales service

For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
 468 rue Jacques Monod
 CS 21900
 27019 EVREUX CEDEX FRANCE

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediately to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.

** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux