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Electromagnetism

**Réf :
292 075**

Français – p 1

English – p 11

Version : 8006

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1 General

1.1 Purpose of the instrument

The reversible demonstration alternator is designed to experimentally illustrate the influence of relative movement of a magnet and a coil to generate a voltage.


The variable voltage in time is obtained by rotating the coil relative to the two magnets or rotating the two magnets relative to the coil.

The reversible demonstration alternator can also be used to show that the maximum value of the voltage and the period depend on the rotation speed.

The alternating nature of the voltage generated may be demonstrated simply using a LED connected to the coil terminals.

A more detailed study can be carried out using an oscilloscope or an acquisition interface connected to a computer.

1.2 Description

1.4 Use

The alternator is positioned in one direction, and the upper plate to which the two magnets are fixed is free to rotate relative to the coil.

The alternator is positioned in the other direction, the upper plate to which the coil is fixed is free to rotate relative to the magnets.

The upper plate is rotated:

- either by pushing one of the four stops fixed to this plate with your finger, and holding the lower plate fixed.
- or by moving this plate uniformly, using a pen located in the orifice provided for this purpose.

1.5 Technical characteristics

The device is made from transparent polycarbonate. The two coaxial plates are mounted on a ball bearing.

The coil inductance is 650 mH and the magnets are made of ferrite.

The coil ends are fixed to two pairs of $\varnothing 4$ mm safety bushings in parallel at a pitch of 38 mm compatible with UME dipoles.

The two pairs of bushings in parallel can be used to connect a dipole (for example a LED) to the coil so that it is always visible regardless of the position of the reversible rotating alternator (see images above).

2 Usage

2.1 Production of a variable voltage

2.1.1 Demonstrating the alternating nature of the voltage produced

The alternator is positioned so as to rotate the magnets.

A LED is connected to the coil terminals.

The magnets are rotated relative to the coil, either by pushing them or using a pen located in the orifice of the upper plate.

If the rotation speed is sufficient, it will be seen that the LED flashes. The LED is special in that it only lights up when current passes through it in one direction, consequently it can be concluded that rotating magnets relative to the coil generates an alternating voltage at the coil terminals that only generates a current during one alternation (the LED lights up).

If the LED blinks when the rotation speed is sufficiently high, this means that the maximum value of the variable voltage depends on this rotation speed.(see section 2.3.2)

2.1.2 Demonstrating the relative movement of magnets and the coil to produce a voltage

The alternator is turned over, and the LED is connected to the bushings of the upper plate.

The coil is rotated relative to the magnets.

The LED blinks when the rotation speed is sufficient.

It can be concluded that an alternating voltage is generated at the coil terminals.

2.2 Power level indicator

The power level indicator connected to the alternator coil terminals demonstrates that the available electrical power* depends on the rotation speed.

* **Note** : The electrical power output by the alternator depends on the receiver connected to its terminals.


2.3 Displaying the voltage generated using an interface

The coil is connected to an input of the chosen interface (Visio, Primo or VTT provided with a voltmeter adapter), to display the voltage generated by the rotating alternator.

2.3.1 Voltage when the alternator is not connected to a LED

The alternator is rotated at no load (without LED), by pushing it with a finger.

The voltage acquisition over a relative long period gives the following result:


Observation of the variable voltage with time demonstrates two important phenomena :

1. The maximum value of the voltage decreases when the rotation speed reduces.
2. The period (or pseudo period) increases when the rotation speed reduces.

2.3.2. Voltage when the alternator is connected to a LED

The alternator connected to a LED is rotated by pushing on it.


Observation of variable voltage with time demonstrates the behavior of the LED:

1. The LED blinks because it only lights up during one alternation. In the experiment made, the LED lights up during the positive alternation (if the LED connection were inverted, it would light up during the negative alternation). Like all diodes, the LED only allows current to pass in one direction.

2. The LED blinks if the value of the voltage at its terminals is sufficiently high, which demonstrates the existence of a «threshold voltage» as for all diodes.
3. The LED limits the voltage in the conducting direction.

2.3.3 Voltage as a function of the relative position of the coil and magnets

When the alternator turns slowly, the plot of the voltage curve as a function of time and the position of the rotor (the rotating part of the alternator) relative to the stator (the fixed part of the alternator) can be observed simultaneously.

It is found that:

1. the maximum and minimum voltages are obtained when the coil axis is perpendicular to the straight line joining the magnets.
2. the voltage is zero when the coil axis is coincident with the straight line joining the magnets.

Note : With the Visio (and Visio plus) console, the alternator can be filmed in movement during acquisition of the voltage. Thus, the variation of the voltage during time and the position of the coil relative to the magnets can be observed, image by image and simultaneously.


Starting from these observations, the alternation of the voltage and the direction of the winding in coil turns relative to a magnet can be related to each other.

The concept of magnetic flux can then be introduced, depending on the study level.

2.3.4 Influence of the coil rotation direction on the voltage produced

The coil is connected to an interface connected to the computer. The polarity of the magnets is identified using a magnetized needle and the south pole of a magnet is placed facing the face of the coil connected to the safety bushings.

Acquisition is started and then the magnets are made to rotate briefly relative to the coil, firstly in one direction and then in the other direction. The results obtained are as follows :

The rotation direction has no influence on the shape of the signal. In both cases, the voltage produced as soon as the magnets start rotating is positive increasing and then decreasing.


2.3.5 Influence of the polarity of the magnets and the direction of winding of the wire that forms the coil

The coil is connected to an interface connected to the computer. The polarity of magnets is identified using a magnetized needle.

In a first experiment, the south pole of a magnet is put into place facing the face of the coil connected to the safety bushings.

In a second experiment, the north pole of a magnet is put into place facing the face of the coil connected to the safety bushings.

For each experiment, the acquisition is started and then the magnets are made to rotate briefly relative to the coil. The results are as follows:

In the first experiment, the south pole (left magnet) « sees » the winding of the coil wire in the anti-clockwise direction.

In the second experiment, the south pole (right magnet) « sees » the winding of the coil wire in the clockwise direction.

It is found that the voltage produced is positive in the first experiment and negative in the second.

In conclusion, when the winding direction of the turns in a coil facing a magnet pole (north or south) changes, the voltage generated is inverted.

3. Theoretical alternator design

3.1 Statement of Lenz's law

The direction of the induced current is such that its effects oppose the cause that created it.

3.2 Magnetic flux

Let N be the number of turns per row on the length L of the coil ($L = 38$ mm).

(The diameter of the copper wire that forms the coil is 0.2mm, while the number of turns per row is $N \approx 190$).

Let i be the number of the row of turns with radius r_i (for $i = 1$, $r_i = 4.5$ mm)

i varies from 1 to n_T (where $n_T \approx 33$)

The expression of the magnetic flux is $\Phi = N\pi \sum_{i=1}^{i=n_T} r_i^2 \vec{B} \cdot \vec{n}$.

If we let $N\pi \sum_{i=1}^{i=n_T} r_i^2 = S_T$, be the sum of the areas of the turns making up the


coil, the expression for the flux becomes $\Phi = S_T \vec{B} \cdot \vec{n}$.

3.3 Arbitrary choice of the normal to the coil

The positive direction is chosen to be the turn winding direction. The right hand rule gives the direction of the normal. (See diagram below)

3.4 Electromotive force (e.m.f) of induction : e

Any variation in the flux generates an e.m.f given by the


Coil seen from the bushings side

expression $e = -\frac{d\Phi}{dt}$.

3.5 Study of the voltage at the coil terminals at no load (electromotive force) during time

Consider the coil in rotation relative to the magnets:

In the following study, we refer to the face of the coil from which the normal originates.

Some positions of the coil as it rotates relative to the magnets are marked at different times. The views are seen from above and the explanations are given with reference to:

1. Lenz's law
2. the magnetic flux and variations of the flux

<p style="text-align: center;">①</p> <p>Lenz's law The coil face moving away from the south pole of the left magnet is a north face, to oppose this separation: $e > 0$</p> <p>Magnetic flux The flux is positive decreasing.</p>	<p style="text-align: center;">②</p> <p>Lenz's law The coil face remains a north face to move away from the left magnet and move towards the right magnet. $e > 0$</p> <p>Magnetic flux The flux is null.</p>
<p style="text-align: center;">③</p> <p>Lenz's law The coil face moving towards the north pole of the right magnet is still a north face and opposes this approach. $e > 0$</p> <p>Magnetic flux The flux is negative decreasing.</p>	<p style="text-align: center;">④</p> <p>Lenz's law The coil face does not have a magnetic face. $e = 0$</p> <p>Magnetic flux The flux is minimum.</p>

⑤	⑥
<p>Lenz's law The coil face moving away from the north pole of the right magnet is a south face, to oppose this separation. $e < 0$</p> <p>Magnetic flux The flux is negative increasing.</p>	<p>Lenz's law The coil face remains a south face to move away from the right magnet and move towards the left magnet. $e < 0$</p> <p>Magnetic flux The flux is null.</p>
⑦	⑧
<p>Lenz's law The coil face moving towards the south pole of the left magnet is a south face to oppose this approach. $e < 0$</p> <p>Magnetic flux The flux is positive increasing.</p>	<p>Lenz's law The coil face does not have a magnetic face. $e = 0$.</p> <p>Magnetic flux The flux is maximum.</p>


The variations of the induction e.m.f. can be summarized in the following table:

positions	①	②	③	④	⑤	⑥	⑦	⑧
Flux Φ		0		mini		0		maxi
$e = -\frac{d\Phi}{dt}$	> 0	> 0	> 0	0	< 0	< 0	< 0	0

3.6 Experimental verification

The experiment is carried out under conditions described and illustrated above, after marking the polarities of magnets using a magnetized needle. The coil is connected to an interface connected to the computer ($u_a = e$). The coil is rotated using a pen located in the orifice provided for it, in the upper plate, at the same time as the acquisition is started.

The experimental verification gives:

Note : Inflections in the curve at the maximum and minimum points are essentially due to the dimensions of the coil. Furthermore, the magnetic field in which the coil is immersed is not strictly uniform.

4 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service *

+ 33 (0)2 32 29 40 50**

* Call cost not included

** Only for call from foreign countries