Appareil de chauffage

Heating drying oven

**Ref :
535 015**

Français – p 1


English – p 8

Version : 8006

Etuve 20 litres

Oven 20 litres

1 Schéma

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
FRANÇAIS

1 Diagram

1. Case
2. Door
3. Door Knob
4. Panel
5. Temperature Controller
6. Exhaust flap Adjustmen
7. Power Switch of Indicator

2 Scope of application

The oven can be used for such tests as drying, curing, melting, sterilizing in industrial and mining enterprises as well as in laboratory and scientific research institutes.

3 Technical features

Model	DHG-9023A
Mains voltage	220 V, 50 Hz
Temperature range	Ambient temperature +10-200°C
Temperature fluctuation	1°C
Heater power	530 W
Dimensions of the working chamber (mm)	300 x 280 x 275

4 Structure

1. The oven case is made of rigid welded steel sheet. The working chamber is made of stainless steel. The heat insulation material adopts fibre glass which is in the case jacket. The case door is made of double glass. The entire unit consists of three parts-hot air circulation system, of an electric heater and a temperature controller system.

2. The hot air-circulation system is composed of single-phase motor, centrifugal vane wheel and wind direction. The electric heater is installed directly in the chamber. The fan motor keeps constant air circulation. It can pick up speed of rising temperature and also keep air temperature fluctuation constant in the chamber.


3. A high-precision temperature controller is installed in the oven. The temperature controller is equipped with a digital display for the setting value and actual value of temperature.

4. In case the temperature in the chamber is 10°C higher than the preset temperature, the temperature controller will switch on the alarm. At the same time, the power supply of heater will be cut off.

5. The air in the chamber can be adjusted with the air flap.


6. The microprocessor intelligent temperature controller is equipped with self-adjust PID method. Compared to normal PID controlling method, self-adjust PID method is of higher accuracy and rapid controlling, the setting value of the temperature and the actual temperature are all digitally displayed. The controller can be easily to operate though a soft-touch pad. One of the standard functions is the alarm that is switched on when temperature exceeds the upper temperature limit.

5 The function panel of the controller and the operation procedure

1. Actual temperature display (PV)
2. Heating indicator light
3. Upper limit alarm indicator light
4. Self-adjust indicator light
5. Setting key
6. Shifting key
7. Decreasing key
8. Increasing key
9. Temperature setting display (SV)

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Remark : Except the setting of temperature, all the other data of setting have been locked. If you wish to revise the new data. First cancel the command for locking out modification of the default parameters.

6 Operation

1. Put the sample into the chamber. Close the door and turn the 'Exhaust flap' to 'Z' position.

2. Turn the power switch to « 1 ». The indicator of the power supplier will turn on and the temperature will be digitally displayed on the controller.

3. Temperature setting

After having switched on the unit, if the set temperature is the same as the needed working temperature, then it is not necessary to set it again. Otherwise, it is necessary to reset the working temperature. First, press the function key 'SET' on the control panel to enter the temperature setting program, till the setting temperature display SV is flashing. Then press the 'SHIFTING KEY' and the 'INCREASING KEY' or the 'DECREASING KEY' to set the needed working temperature. Finally, press again the 'SET' key for 5 seconds to confirm the setting parameters.

For example, the needed working temperature is 150.0°C, and the former setting temperature is only 86.5°C. First, press the function key 'SET' on the control panel to enter the temperature setting program, till the setting

temperature display SV is flashing. Then, press the 'SHIFTING KEY' to move the cursor to the third digital position before the base point, then press the 'INCREASING KEY' to change the display from « 0 » to « 1 ». Then press again with the « 'SHIFTING KEY' » and the « INCREASING KEY » or the « DECREASING KEY » to set each digital position until 150.0°C appears on the display as needed working temperature. Finally, press again the 'SET' key for 5 seconds to confirm the setting parameters.

All the following functions are pre-set by default. They can only be modified after cancelling locking out of parameter modification = "LCK" function.

4. Upper limit temperature alarm setting

Before delivery, the upper limit temperature alarm setting of all units was set at 10 °C .

Therefore, normally it is not necessary to re-set it again. If it is necessary to re-set it again, then first press the function key « SET » for 5 seconds, the controller enters the upper limit temperature setting program AL1. Press the 'SHIFTING KEY' and the 'INCREASING KEY' or the 'DECREASING KEY' to set the demanded upper limit temperature alarm setting. Finally, press again the 'SET' key for 5 seconds to confirm the setting parameters. (See 6.10)

5. Temperature display correction

Before delivery, all units were carefully tested, normally it is not necessary to correct the temperature display. If the working environment is not good (ie if the environment temperature is too high or too low), the temperature displayed could be different from the actual temperature, then it is necessary to correct it. These are the steps to follow : first press the function key 'SET' for 5 seconds to enter the parameter setting cycle AL1. Then keep on pressing the 'SET' key, till the PV display 'SC' as temperature display correction. Then press the 'SHIFTING KEY' and the 'INCREASING KEY' or the 'DECREASING KEY' to set actual temperature display. Finally, press again the 'SET' key for 5 seconds to confirm the setting parameters.

6. After setting, all parameters will be stored for a long period of time. The unit will start heating up, the heating indicator light will turn on. While drawing closer to the setting temperature, the heating indicator light will keep on flashing, till the unit will reach the setting temperature.

7. The selection of heating time is determined by the degree of humidity of the samples. If the sample has a high degree of humidity, turn the « Exhaust flap » to ≡ position to ventilate the wet air out of the unit.

8. After you have completed the whole heating procedure, close the ' Exhaust flap' if the samples were not take out at once. Otherwise, still keep the 'Exhaust flap' open, and turn the power switch to ' 0 ' position. Be careful for the hot samples while taking them out.


9. If more accuracy is needed in temperature control and deviation, then the self-adjust PID function can be useful. When the temperature of the units is almost close the to setting temperature, first press the function key 'SET' for 5 seconds to enter the parameter setting cycle ' AL1' . Then keep on pressing the 'SET' key, till the PV display 'ATU' . The SV display '0000' ,then press the 'INCREASING KEY' to make SV display as '0001' .Finally press again the 'SET' key for 5 seconds to confirm the setting parameters. Then the self-adjust indicator light will turn on. The unit will enter the self-adjust PID control.

10. Unlocking to modify the pre-settings :
 Press the "SET" button for 5 seconds, the "AL1" menu will be displayed, press "SET" by scrolling the various functions until "LCK". Using the buttons "Δ" or "∇" displayed, correct the value to "0". Press "SET" to confirm. Locking out is now cancelled. The functions can be modified. To relock the commands and so avoid any accidental changes ; repeat the above procedure and select a value other than "0" for the "LCK" function.


7 Please pay attention that

1. The outer casing of the heater must be effectively grounded so as to ensure safety in use.
2. In case any alarm is switched on, please switch off power at once, and just trace out the cause, and do not use the device until it is repaired.
3. Don't place any combustible and explosive goods in chamber or near the oven.
4. If the unit will not be used for a long time, it should be covered with a plastic film dust cover and be placed it in a dry and well-ventilated room.

8 Block-diagram of the electrical principle

9 Wiring diagram

10 Trouble shooting


Failure	Cause	Shooting method
1. Display failure	<ol style="list-style-type: none">1. The device hasn't been plugged in or the electric cable or plug is out of work.2. Fuse is blown.3. The power switch is not in the « on » position.	<ol style="list-style-type: none">1. Examine & repair the plug.2. Replace the fuse.3. Place the switch in its right position.
2. Temperature will not rise	<ol style="list-style-type: none">1. The preset temperature is too low.2. The heater is out of order.	<ol style="list-style-type: none">1. Adjust the preset temperature.2. Replace the heater.
3. Overtemperature alarm	<ol style="list-style-type: none">1. Operation is wrong.	<ol style="list-style-type: none">1. First switch off power and open the door for a moment, then restart the machine.

11 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
0825 563 563

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service *

+ 33 (0)2 32 29 40 50**

* Call cost not included

** Only for call from foreign countries