# **Electricité**

Electrostatique

# ***Electricity***

*Electrostatics*

**Ref :  
272 006**


**Français – p 1**

**English – p 5**

**Version : 9003**

**Machine de Wimshurst**

***Wimshurst machine***

## 1 Principe – description

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
**NOTES**

## 1 Principle – description


The Wimshurst machine consists of 2 isolating plates of identical size (30 cm) D<sub>1</sub> and D<sub>2</sub> separated by a few millimetres. They rotate in opposite directions, and constitute the charge carriers.

There is a transverse, diametric conductor equipped with small metal brushes rubbing against the tin strips in front of each disk.

Fixed metal rakes, P and P' communicate with the capacitor plates, and the terminals of the machine.

When one of the strips is charged + on D<sub>1</sub> for example, the opposite strip on disk D<sub>2</sub> gets charged negatively. By influence, the disks get charged according to the signs indicated in the diagram hereof.

These charges driven by the rotation maintain the distribution as shown on the next page.

The positive charges are captured by the collector P' and the negative charges by collector P. They are then channelled towards the 2 conductors with spherical tips at the end that act as spark gaps.

By introducing the Leyden jars in the circuit, with the connections of the device we increase the electrical energy discharged during sparking.

**NB:** we have supposed that the machine is already operating. In fact the machine self-starts, either because there are potential differences between the various parts, or because the friction of brushes is enough to create them. (Start with a charged ebonite stick, otherwise.)

## 2 Experiments

### 2.1 Check the position of various brushes

- The brushes fixed on the two diametric conductors must rub lightly against the metal strips.
- The “rake” brushes must be close to the strips. It is not necessary for them to rub against the strips for the charges to be transported.
- The contact brush → collector distance will be adjusted so that the charge carried by each disk is maximum (see diagram).

The influence  $D_1 \leftrightarrow D_2$  will then be maximum. The brushes will thus be close to the collectors at a distance greater than the starting voltage of the spark gap.

### 2.2 Open circuit (the two connections near the collectors raised)

- Rotate the handle clockwise, with the spark gaps 2 to 4 cm from one another. Crackling is heard due to discharges of very low intensity between the spark gaps. These discharges are not very visible. Ozone is emitted.

### **2.3 Repeat the experiment by closing the circuit on the Leyden jars**

- Rotate the handle forcefully.
- When the two jars are charged, a clearly visible spark appears, particularly if the connecting strip at the base of the device closes the circuit.


## **3 After-Sales Service**

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

**JEULIN - TECHNICAL SUPPORT**  
**Rue Jacques Monod**  
**BP 1900**  
**27 019 EVREUX CEDEX FRANCE**  
**+33 (0)2 32 29 40 50**

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

## Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

### Free service \*

**+ 33 (0)2 32 29 40 50\*\***

\* Call cost not included

\*\* Only for call from foreign countries