Didactique

Kit de chimie

Teaching manual

Chemistry kit

Ref :
106 057

Français – p 1

English – p 5

Version : 1005

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Student chemical garden kit

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Various salt crystals are introduced in a container that contains sodium silicate. Coloured membranes grow from the surface of each crystal. The more hygroscopic the crystal, the faster the membranes grow.

1 Composition of the demonstration set

- 5 bottles that can contain 50 ml of dilute sodium silicate Na_2SiO_3 solution.
- 1 bottle containing 250 ml of dilute sodium silicate Na_2SiO_3 solution.
- A set of 15 bottles containing the following salts: $\text{CuSO}_4, 5\text{H}_2\text{O}$; $\text{FeCl}_3, 6\text{H}_2\text{O}$; $\text{FeSO}_4, \text{H}_2\text{O}$.

1.1 Precautions

- Sodium silicate Na_2SiO_3 is a strong base.
- $\text{CuSO}_4, 5\text{H}_2\text{O}$ CAS: [7758-99-8] Xn R: 22-36/38-50/53
S: 2-22-60-61
- $\text{FeCl}_3, 6\text{H}_2\text{O}$ CAS: [10025-77-1] C Xn R: 22-34
S: 2-7/8-28
- $\text{FeSO}_4, 7\text{H}_2\text{O}$ CAS: [7782-63-0] Xn R: 22-41
S: 2-26

1.2 Waste disposal

Allow the solution to decant or filter. Eliminate with excess water or neutralise with hydrochloric acid until the solution clears.

Place solid waste in an appropriate recycling container.

2 Operating procedure Presentation

Shake the bottle to obtain a homogeneous solution.

Fill each bottle with 50 ml of solution, place the individual bottle on a flat surface and wait for the solution to become stable.

Open the bottle and introduce the salts one by one sprinkling them over the bottom surface of the bottle. The colloids have a specific colour and size depending on the salt deposited. The higher the particle size of the salt, the faster it reacts.

Close the bottle once the crystals are inside and do not move the chemical garden. Movement may break the colloids formed.

The progression of the salts may continue for several days.

3 Discussion


The formation of colloid membranes of the chemical garden is possible due to two combined phenomena: and buoyancy.

The presence of sodium silicate allows the formation of metal silicates:

When they form, the metal silicates have a semi-permeable barrier that allows obtaining coloured membranes.

3.1 Explanation of the first phenomenon


When the semi-permeable barrier has formed (1), water comes in through osmosis but cannot leave (2). Therefore, the salt concentration decreases and the water volume increases. The pressure increases until the barrier bursts (3). When it bursts, an amount of salt with a concentration lower than the initial concentration C_1 is expelled. These salts react with the silicates again to form a new barrier and so on (4 and 5)

Why does it progress towards the top?

3.2 Explanation of the second phenomenon

When the barrier forms, the density inside it is higher than the density of the medium (1). As the water enters by osmosis, the salt density decreases because it dissolves (2). When the membrane breaks up (3), the salt density is lower than that of the medium. Due to the buoyancy principle, the low density solution goes up to the surface explaining the F force applied from bottom to top.

4 After-Sales Service

For any repairs, adjustments or spare parts, please contact:

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediately to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service *

+ 33 (0)2 32 29 40 50**

* Call cost not included

** Only for call from foreign countries