**Méthode ESAO<sup>®</sup>**

***ESAO<sup>®</sup> method***

**Ref :  
472 067**

**Capteur Electrophy**

***Electrophy sensor***

Français – p 1

English – p 9

Version : 8012

# 1 . Généralités

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
# 1 General

## 1.1 Purpose of the instrument

The Electrophy sensor is designed to make data acquisitions in human electrophysiology. It performs two functions and can

- in the Electrophy position, perform raw signal acquisitions (ECG, myotatic reflex, electromyography)
- in the heart rate meter position, make a heart rate measurement.

The function is chosen using the cursor on the top of the sensor

It is composed of a sensor to connect to your interface (PRIMO, VTT, SYNEO)

- a 3-wire cable with 3 colour plugs to connect the electrodes.
- 3 wrist electrodes for monitoring the change in the heart rate.

The product is delivered with a flask of gel to enable better measurement.

## 1.2 Recognition by the PRIMO console

With the PRIMO console, the Electrophy sensor is recognised by the generalist scientific studio software starting from version 2.2. and by dedicated scientific studio software starting from version 2.3.

### 1.3 Recognition by the SYNEO console

Your SYNEO console natively recognises the Electrophy sensor.

### 1.4 Recognition by the VTT console

Recognition by the VTT console requires your system to be updated. You can do this update from the [www.jeulin.fr](http://www.jeulin.fr) internet site, or using the disk delivered with the sensor.

This update is made when your console is connected to the PC using the Espace-VTT software:

- Connect the console to your PC
- Open the VTT space by clicking on the  icon on your workstation
- Click on install 
- Install the Biosflash software.

**You can then use it with VTT:**

- in standalone mode
- or
- connected to a computer (a scientific studio acquisition software is then necessary).

## 2 Example applications

### 2.1 Performing an ECG


#### 2.1.1 Placement of electrodes

When used in secondary schools, the easiest way is to ask pupils to place their electrodes on the inside of the wrists.

##### Placement on the wrist:

- Apply a little conducting gel on the wrist.
- Put the straps into position and tighten them so as to obtain a thin coat of gel between the metallic part and the skin.
- The yellow electrode is placed on the left side,
- The red electrode is placed on the right side,
- The violet electrode is placed on the right ankle.

It is recommended that the wire should be looped under the strap to prevent the connector from moving during displacements (see photo on next page).

In this case, when taking a measurement with electrodes on the wrist, it is recommended that the patient should be resting, at least for the upper part of the body, otherwise the system will measure the muscular electrical activity in the arm.


### Placement on the thorax:

Placing electrodes on the thorax gives better quality results. Furthermore, in this case, the signal shape is genuinely the shape of an ECG with PQRS waves.


- The yellow electrode is placed on the left side.
- The red electrode is placed on the right side.
- The violet electrode is placed on the back.


## 2.1.2 Software parameter settings

- Drag the icon representing the Electrophy sensor (in the signal position), into the ordinates
- Drag the alarm clock icon representing the time, into the abscissa.

- Set time parameters as follows:


For a better representation, it is recommended that points should be connected. Proceed as follows for this operation:

- Click on the Electrophy icon and then on the display tab and then on Connect

- Start acquisition by clicking on the green light

## 2.1.3 Results

## 2.2 The Myotatic reflex

### 2.2.1 Placement of the sensor

Sticky electrodes are necessary for the myotatic reflex.

Arrange them in a triangle at the top of the calf, with the violet electrode towards the bottom.

The pupil must keep still with his knee on a chair to minimise electrical activity of muscles necessary to stay in position.


Connect the Martoreflex onto the synchro input of the interface (red and black jacks).

### 2.2.2 Software parameter settings

- Drag the icon representing the Electrophy sensor into the ordinates
- Drag the alarm clock icon representing the time into the abscissa
- Set time parameters as indicated on capture 1 on the next page

The beginning of the data acquisition has to be synchronised for the myotatic reflex.


This is done by clicking on the synchronisation tab and choosing Synchro input in the pull down menu ❶ and select decreasing ❷ as shown on the capture 2 below

- Start acquisition by clicking on the green light 

A message is displayed to tell you that the machine is waiting for synchro. Use the previously connected Martoreflex and tap the foot tendon with it.

## 2.2.3 Results

## 2.3 Monitoring the heart rate

### 2.3.1 Placement of electrodes and hardware parameter settings


- Put the electrodes into place as for the ECG (See chapter 2.11.)
- Put the sensor into the rate position.
- Adjust the gain using the potentiometer until the red LED corresponding to the heart beats flashes on the sensor.

If the LED is off, the gain needs to be increased and if the LED is continuously on, then the gain must be reduced.  
There is a short latency time between when the gain is modified and its impact on the behaviour of the LED


### 2.3.2 Software parameter settings

Drag the icon representing the Electrophy sensor (in the signal position), into the ordinates


- Drag the alarm clock icon representing the time, into the abscissa.

- Set the time parameters as follows:

### 2.3.3 Results

## 3 Technical characteristics

### 3.1 Usage precautions and standards

- The product has a CF type insulation (degree of protection against electrical shocks controlling the product insulation).

5335

- The product must be used only with Jeulin consoles and their mains adaptors.
- Operation guaranteed only with bracelets supplied by Jeulin,
- Conducting gel must be used
- The computers used must comply with standard IEC950.

### 3.2 Usage conditions:

Working temperature:	10 to 50°C
Pressure:	from 700 to 1300 hPa
Relative humidity:	30% to 95%

### 3.3 Electrode characteristics

Anallergic 302 stainless steel  
Anallergic conducting bracelet

## 4 After-Sales Service

For any repairs, adjustments or spare parts, please contact:

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

## Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

### Free service \*

**+ 33 (0)2 32 29 40 50\*\***

\* Call cost not included

\*\* Only for call from foreign countries