Géologie

Tectonique

Geology

Tectonics


**Réf :
507 037**

Français – p 1

English – p 7

Version : 7005

Machine tectonique élève
Student's tectonic machine

1 Principe

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
NOTES


1 Principle

A transparent tank is divided by a wall moved by a piston. This device allows reproducing through the pressure exercised on materials (sand, plaster etc.) in super-imposed layers, the **main consequences in distortions linked to tectonic movements**, notably phenomena of compression, folding, faults, reverse faults.

2 Composition

The tectonic machine for students is composed of:

- A tank with a cover, in polypropylene, dimensions 140 x 100 x 70 mm
- A pushing wall
- A piston
- A tamper to compact the materials

3 Assembly

Remove the cover and other components from the pierced tank.

Pass the end of the piston through the orifice of the tank and push into the housing of the pushing wall which has previously been inserted in the tank. Position the wall carefully, with the rounded corners at the bottom and the sharp corners at the top.

4 Use

4.1 Materials and preparation

The principle consists in distorting layers of materials of different nature and/or colour:

- White plaster
- Coloured plaster
- Natural sand (natural or artificial such as micronised plastic), of different colour and grain size.
- Other materials, such as washing powder

The different materials can be alternated.

The interposition of a fine sandy layer (plastic or washing powder) between two layers of plaster can act as a non-adhesive layer (soapy layer).

Materials can be recovered at the end of the operation for mixing and then used as basic materials for the next operation.

Use of plaster or cement allows retaining the models constructed. After drying out they can be cut into cross-sections.

Insertion of layers of materials:

Increase the working distance by withdrawing the pushing device on the piston side.

The various layers of plaster or other materials are laid regularly, each layer being of uniform thickness throughout. Compact the material using the tamper, in particular around the edges in contact with the walls of the tank to create clear lines.

To optimise tamping, correctly orient the tamper by placing the rounded corners towards the end of the tank and the sharp corners next to the pushing device.

Irrespective of the layout of materials, each compression operation will generate different figures depending on the nature, distribution and thickness of the layer.

There is no point in the upper layer exceeding mid-height of the tank since compression will lift up the material which could then spill out.

As a general rule, fine grain size (plaster) and strong compression promote the obtaining of faults. Larger grains which do not stick together promote the formation of folds.

4.2 Operation

With the layers of material in place, in one hand hold the end of the tank containing the layers to be compressed.

With the other hand, grip the end of the piston and push it against the palm of the hand.

Apply very gradual pressure to cause slow movement of the pushing wall. Compression of the structures begins. Observe the behaviour of different layers which are shortened, causing them to rise up and become distorted, whether in the form of breaks or folds, depending on the type of materials used.

During these movements, it is possible to stop progress at several stages to observe the distortions and then carry on to the required stage.

To preserve the model, slowly immerse the tank in a recipient full of water, leaving the water to enter slowly through the piston orifice until the materials are soaked. Slowly remove the tank and leave it to set before removal of the materials.

4.3 Results

The structures which can be observed are as follows:

- Folds
- Faults
- Tectonic flaked faulted folds
- Overlaps
- Inversion of series and contact between terrains of very different ages
- Reverse drift
- Elongation fractures

5 Maintenance

The tank is made of polypropylene which has the advantage of greater scratch resistance than many other plastics. However, use of mineral materials, notably sand, does not totally allow avoiding scratching. It can be limited or even avoided by using materials less abrasive than sand, such as plaster, washing powder, plastics...

Polypropylene also provides good chemical resistance to a number of organic liquids. Nevertheless avoid contact with organic solvents. Wash with water or dishwashing products. Machine washing is possible.

6 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediately to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service * :
+ 33 (0)2 32 29 40 50

** Call cost not included*