# **Didactique**

Thermochimie

# ***Teaching Presentation***

*Thermo-chemistry*

Réf :  
253 059

Français – p 1

English – p 7

Version : 8005

**Etoile 4 métaux Al, Cu, Fe, Zn**

**4 metal star Al, Cu, Fe, Zn**

## 1 Description

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**NOTES**

## 1 Description

The 4 metal star (Al, Cu, Fe, Zn) provides a simple, clear method for introducing students to the thermal conductivity of metals. The apparatus allows comparing the heat conduction properties of 4 different metals in a single experiment, that is aluminium, copper, iron and zinc.

The 4 metal star is constituted of 4 blades (Al, Cu, Fe, Zn) of exactly identical dimensions and thickness, attached to a circular plate which acts as a heat distributor. The results obtained are thus scientifically rigorous.

### 1.1 General

#### 1.1.1 Description


Photo 1

- |                            |                 |
|----------------------------|-----------------|
| 1 – Metal star handle      | 3 – Metal blade |
| 2 – Heat distribution disk | 4 – Crucible.   |

Table 2

#### 1.1.2 Additional equipment

Use of the 4 metal star requires a minimum amount of equipment.

- A support,
- A tea light,
- 4 pieces of paraffin wax (taken from the tea light),
- A cigarette lighter or box of matches.

The 4 metal star can be used more safely with its special support (*Photo 3*) which is sold separately.

4 metal star support P/N 253 065


Photo 3

- |                     |                            |
|---------------------|----------------------------|
| 1 – Retaining slide | 2 – Housing for tea light. |
|---------------------|----------------------------|

## 2 Precautions for use

### 2.1 Warning

JEULIN cannot be held liable in the case of an accident which occurs on account of non compliance with these safety instructions.

### 2.2 Risks of burns

During the experiment the 4 metal star is heated to a high temperature and care must be taken to prevent any risk of burns which may occur in the case of direct contact with the apparatus.

The 4 metal star must be handled using pincers or heat-proof gloves both during and after the experiment.

If necessary after the experiment the 4 metal star can be cooled by immersing it in a bath of cold water.

Those using the 4 metal star must wear a cotton overall and tie back their hair.

### 2.3 Risk of fire

The 4 metal star is designed to be heated with a tea light type candle. Do not use gas or liquid fuel heaters or burners for the operations.

Before carrying out the experiments, ensure no flammable or combustible products and no fuel are close to the work bench.

Always work on a non-combustible, clean, empty surface. Avoid sheets of paper which may be on the work bench. If applicable remove these away from the apparatus.

Walkways around the work bench must be clear so that people can circulate freely.

## 3 Operating method

### 3.1 Setting up the apparatus

Attach the 4 metal star at a height of 4 cm above the bench surface. To do this a stable support with a weight and a clip can be used. Use of the optional metal star support is nevertheless advised since it ensures satisfactory stability of the assembly and guarantees the 4 metal star is at the correct height.

### 3.2 Assembly on the 4 metal star support

To install the 4 metal star on the support, proceed as indicated below (*diagrams 4 and 5*).


Diagram 4


Diagram 5

### 3.3 Preparing for the experiment

Place the tea light under the 4 metal star centring it in the housing provided (*Photo 3-2*). Centre the middle of the heat distribution disk (*Photo 1-2*) vertical to the centre of the tea light.

Place small pieces of paraffin wax of the same thickness (2 to 3 mm in diameter) in the crucibles located at the end of the metal blades.


### 3.4 Performing the experiment

The experiment consists of comparing the time needed for the paraffin wax to melt in each crucible. Since the paraffin wax has a specific melting temperature, its change of state indicates the time the ends of the metal blades, close to the crucible, reach this melting temperature.

Since the metal blades are of the same dimension, it is possible to carry out a direct comparison of the heat conduction properties of each metal.

The experiment may be completed by measuring the time necessary for the morsels of paraffin wax to melt using a chronometer.

Light the tea light with a cigarette lighter or matches. If the experiment is accompanied by time measurement, start the chronometer that moment.

Morsels of paraffin

Photo 6

Observe the different blades and note the order of melting of the paraffin wax in the various crucibles. If applicable, record the time which has elapsed since the candle was lit.

The metals constituting the metal blades are given below (Diagram 7).


Diagram 7

The thermal conductivity of the various metals is given in table C below. (table 8).

Metal	Thermal conductivity <sup>(1)</sup> W.cm <sup>-1</sup> .K <sup>-1</sup>
<b>Aluminium</b>	<b>2.37</b>
Silver	4.29
<b>Copper</b>	<b>4.01</b>
<b>Iron</b>	<b>0.80</b>
Gold	3.18
<b>Zinc</b>	<b>1.16</b>

(1) to 298,2 K – Source: Handbook of chemistry and physics 61<sup>th</sup> edition  
 CRC Press

Table 8

## 4 Maintenance and storage

### 4.1 Oxidation

The metals constituting the metal blades are raw and untreated to prevent distorting the experimental results. For this reason, they are liable to spontaneous oxidation. This is particularly applicable to the iron blade which may rust during storage.

**Important:**

**On delivery the star may show traces of oxidation. This does not compromise the product quality in any way.**

To restore the initial appearance of the metals, gently scrub the blades using a household pan scour or a steel wool pad.

## 4.2 Cleaning the blades

At the end of the experiment and after cooling the 4 metal star it is necessary to remove deposits of solidified paraffin wax from the crucibles.

**Do not scratch the blades or the crucibles with a metal object or tool. This could cause deep scratches in the metal and permanently damage the apparatus.**

To remove the paraffin wax immerse the 4 metal star in a bath of very hot water and when the paraffin wax softens, remove by rubbing the blades and crucibles with absorbent kitchen paper.

Dry the 4 metal star carefully with a dry cloth to avoid oxidation of the metals.

## 4.3 Storage

Store the 4 metal star away from sources of damp. If during storage traces of oxidation appear, refer to paragraph 4.1.

## 5 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

**JEULIN - TECHNICAL SUPPORT**  
Rue Jacques Monod  
BP 1900  
27 019 EVREUX CEDEX FRANCE  
0825 563 563

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

## Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

**Free service \* :**  
**+ 33 (0)2 32 29 40 50**

*\* Call cost not included*