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Conductimétrie

Chemical measurement

Conductimetry

Ref :
701 329

Français – p 1

English – p 5

Version : 7006

Sonde conductimétrique CTM
CTM conductimeter probe

1 Description

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FRANÇAIS

1 Description

The CTM conductimeter probe equipped with graphite electrodes is robust and easy to maintain.

It can be used:

- Connected to the conductimeter CTM Initio[®] Ref. 701 327.
- Mounted with traditional measurement instruments (ammeter and voltmeter) and a sinusoidal voltage generator (500 Hz min.).


Diagram 1


2 Technical characteristics

Electrodes	Graphite
Cell constant	Theoretical ⁽¹⁾ : 0,393 cm Experimental ⁽²⁾ : Measured at 500Hz : 0,69 - 0,75 cm Measured at 5kHz : 0,64 - 0,69 cm
Minimum height of liquid required	20 mm
Dimensions (mm)	Ø body: 12; Ø tip: 16 Length: 140
Connection	Banana safety plugs (Ø 4 mm) with rear outlet
Measurement range: 0 - 20 mS.cm⁻¹	Connected to the CTM Initio [®] conductimeter (working frequency: 5 kHz).
0 - 8 mS.cm⁻¹	In a traditional circuit supplied with sinusoidal voltage of 500 Hz

⁽¹⁾ - Calculated on the basis of relation $K = S/L$. S is the surface of an electrode and L is the distance between the electrodes.

⁽²⁾ - Average values measured on a batch of 10 probes at 500 Hz and 5 KHz.

2.1 Geometry of the probe

2.2 Resistance to aggressive chemical products

This probe is designed for resisting the immersion in saline, acid, basic, diluted electrolytic solutions.

Do not dip this probe in oxidising solutions (potassium permanganate, sodium hypochloride, etc.) or in organic solvents.

For an optimum service life, it is recommended to always wash the probe with distilled water after each usage.

3 Connection

This probe is equipped with a cord provided with banana safety plugs ($\varnothing 4 \text{ mm}$) with rear outlet.

It is connected to the conductimeter CTM Ref.: 701 327 (refer to its manual).

It can also be used with traditional measurement instruments (ammeter, voltmeter) supplied by a sinusoidal voltage generator of 500 Hz min.

4 Stability of measurements

If the probe is not used for a long time and/or if it is stored in a dry ambience, it is recommended, before a series of measurements, to humidify the surfaces of graphite by dipping the electrodes in distilled water for a few minutes (5 to 30 minutes depending on the degree of drying). This ensures the stability of measurements obtained.

5 Maintenance


Rinse the tip of the probe with demineralised water after each usage (do not use the organic solvents), and store in a dry place.

6 After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.

For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediately to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service * :
+ 33 (0)2 32 29 40 50

** Call cost not included*