Optique

Optics

**Réf :
201 018**


Français – p 1

English – p 5

Version : 6101

Diode laser modulable - Vert
Green modulatable laser diode

1 Description

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
1 Description

This 1 mW modulatable laser diode was especially designed for studying information transmission by light, and it complies with the safety regulations established in the National Ministry of Education circular of October 12, 1999. It is equipped with an analog port provided by safety sockets for modulating the light intensity via a TTL signal having a frequency ranging from 0 to twenty kHz.

The apparatus comes with a 5 meter-long fiber optic and a perforated silicon stopper, which allows you to attach the fiber optic quickly and efficiently to the end of the diode.

The diode is also equipped with a $\pm 3^\circ$ vertical adjustment system for the beam, which facilitates the optical alignment for standard uses.

An off/on switch saves the batteries when the diode is not in use.

Attention! Power consumption implies 4h autonomy.

2 How a laser diode works

The laser diode is an optoelectronic device that functions with semi-conductors.

It emits monochromatic, coherent light (an optical power) for carrying a signal containing information. This light has other uses as well.

As with all lasers, the laser diode functions by means of an amplification medium, a resonating structure, and a pumping mechanism.


The advantages of the laser diode are its compactness and the ease of modulation.

3 Technical specifications

Wavelength:	532 nm
Power:	1 mW (class II) Complies with the EN 60825-1 standard
Ø of the beam (at 5 m):	8 mm
Vertical adjustment:	± 3°
Modulation:	TTL voltage (5 V), input ports provided by two 4 mm Ø sockets
Protection:	Against voltage surges of the modulation signal
Power source:	Two AA size 1.5 V (LR6) batteries
Dimensions:	100 x 100 x 40 mm
Weight:	220 g
Running time (battery life):	4 hours of continuous functioning
Fiber optic:	5 meters

4 Functioning

4.1 Power supply

The modulatable laser diode requires two 1.5 V type LR6 batteries.

To open the battery compartment, insert a small screwdriver in the slot on the side of the cover, and pry it open as shown in the illustration on the left.

The green modulatable laser diode requires 2.4V to emit laser beam. For voltage under 2.4V, the fan in charge of cooling box runs even if there's no light emitted.

4.2 Modulation

To modulate the laser beam, send a TTL signal (delivered by a simple function generator) between the 2 black and yellow banana safety sockets located on the side of the apparatus.

The apparatus is electronically protected against voltage surges from the modulation signal.

Frequency should be comprised between 0 and 1kHz to avoid timelag.

4.3 Use with a fiber optic

The laser diode comes with 5 meters of fiber optic for studying the transmission of information by light.

To connect the fiber optic to the end of the laser diode, feed one of its ends through the silicon stopper until several millimeters extend from the other side. Then insert the stopper in the hole provided for this purpose in the laser diode case, making sure that the end of the fiber optic reaches the laser diode.

Turn the unit on and make sure that the green light reaches the other end of the fiber optic.

4.4 Standard use of the laser diode

If there is no modulation voltage applied to the analog input ports, the laser diode functions as a continuous laser source, and can thus be used in standard optical assemblies.

With the $\pm 3^\circ$ vertical adjustment system for the beam, a perfect optical alignment can be achieved with the optical instruments of the assembly.

5 Experiment: fiber optic transmission of a periodic electric signal

In the assembly on the following page, a TTL electric signal is sent to the analog port of the modulatable laser diode, where this electric signal is converted to a light signal. The resulting laser signal is conveyed via a fiber optic to the receiver module of the Fibroptonic box¹ (item no. 202 012 – not included). The resulting electric signal can be observed on the oscilloscope and compared with the original TTL signal.

A phase difference of essentially zero is observed between the incoming signal and the outgoing signal, thus demonstrating the nearly instantaneous transmission of information (in effect, it takes less than 17 nanoseconds for the light to travel through the 5 meter long fiber optic). A weakening of the signal, however, can be observed. This is due mostly to the efficacy of the fiber optic connections to the laser diode on one end and to the receiver on the other end.

¹ The FIBROPTONIC box consists of an assembly of optoelectronic accessories (emitter, receiver, photoresistance, fiber optic) for studying the emission, transmission, and the reception of signals. It can also be used to measure the frequency of various vibratory phenomena.


6 Accessories (not included):

Item name	Item nos.
GTBF	293 449
Fibroptonic	202 012
EVOLIO F3-12 / 1 A power supply	281 001

To observe signals, please connect the laser diode TTL input and fibroptonic output to a data logger or to an oscilloscope. Then, you'll be able to check time lag between emitted and received signals.

Here is the result obtained for a frequency comprised between 0 and 1kHz. Over 1kHz, laser diode response is too long to compare those signals.

7 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
468, rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE
+33 (0) 2 32 29 40 50

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediately to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.

** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux