# Mesure

Initio<sup>®</sup> 2

# Measurement

Initio<sup>®</sup> 2

Réf :  
211 054

Français – p 1

English – p 7

Version : 7105

**Luxmètre Solarimètre Initio<sup>®</sup> 2**  
**Initio<sup>®</sup> 2 Light meter Solarimeter**

## 1. Description

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**0 825 563 563\***

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
## 1. Description

The Initio® 2 Light meter Solarimeter is a device that measures irradiance and illuminance. It can also be used to record data and then restore it later on a computer.

A WI-FI module (ref 451 470) is also available for viewing and recording measurements on a tablet or a PC equipped with WI-FI.

Equipped with an analogue output, the Initio® 2 Light meter Solarimeter can be used with an ExAO interface equipped with direct inputs via banana jacks or via a Groomy® interface and jack.

## 2. Utilisation

Upon first use, it is advisable to completely charge the Initio® 2 battery. Use the  $\mu$ USB (7) supplied in separate packaging to charge the battery. In general, we recommend using the instrument in complete charge-discharge cycles (wait until it is completely discharged before fully recharging)

### 2.1 Start-up

- Press the start button (1) situated on the left side of the instrument.

The measurement displays directly on the instrument's screen.

- To switch off the Initio® 2, briefly press the start button (1) and confirm by pressing the knurled knob (3).

### 2.2 Browsing

Initio® 2 has three navigation buttons:

#### Selection-confirmation (3):

The central knurled knob (3) is to scroll through the menus of the Initio® 2. Briefly pressing the OK button in the centre of the knurled knob confirms the selection.

#### Back (2):

The "Back" button (2) is used to switch from one type of display to another from the main screen in the following order:

- Measurement
- Table
- Curve

Within a menu it is used to go back a step.

#### Menu (4):

This button displays Initio® 2's main menu.

## 2.3 Type of measurement

The Initio® 2 Light meter Solarimeter enables measuring illuminance in *lux* and irradiance in  $W.m^{-2}$

For the light meter, two sizes are available:

- 0 to 2000 lux (measured indoors)
- 0 to 200 klux (measured outdoors)

For the solarimeter, the measuring range is 0 to 2000  $W.m^{-2}$

#### To switch from one measurement to the next:

- Go to the Measurement display
- Circle the name of the parameter on the top right using the knurled knob
- Briefly press the knurled knob

#### To change the calibre of the light meter:

- Set to lux measurement
- Circle the calibre value to the right of the parameter name with the knurled knob
- Briefly press the knurled knob
- 

## 2.4 Data acquisition

Initio® 2 has an internal memory that enables it to save measurements. Two modes are available:

- Manual acquisition

Acquisition is on a point by point basis. Each measurement is manually recorded by the user (see below)

- Time acquisition

Data acquisition takes place as a function of time, according to a duration entered by the user.

### 2.4.1 Manual acquisition

- Press the Menu button (4)
- Using the knurled knob (3), select **Acquisition** and confirm by a short press on the knurled knob (3).
- Select **Manual**
- To record a measurement, select the **Save** icon (white circle) and confirm by briefly pressing the knurled knob (3).
- To stop the recording select the **Stop** icon (black square) and confirm by pressing the knurled knob (3). Confirm the selection in the dialogue box.

### 2.4.2 Time acquisition

- Press the Menu button (4)
- Using the knurled knob (3), select **Acquisition** and confirm by briefly pressing on the knurled knob (3)
- Select **Time**
- Using the knurled knob (3), select the number defining the duration and press on the knurled knob (3).
- Modify this duration by turning the knurled knob (3) then confirm by pressing on the centre of the knurled knob.
- If the data acquisition is to loop (return of the acquisition to zero at the end of a defined period) select **Loop**
- Confirm the final choice by selecting **OK**

Data acquisitions based on time last between 60 seconds and 50 hours. The interval between two measurements adjusts automatically.

### 2.4.3 Recording the data

Once the acquisition is complete (time or manual), it must be saved to be able to retrieve the points on a computer.

- Press the Menu button (4)
- Select **File** and confirm
- Select **Save**

Only a single data acquisition can be saved. By selecting **Open**, you open the last saved data acquisition.

### 2.4.4 Transferring the data to the computer

Once the acquisition is complete, it is possible to transfer the data to a computer to process it:

- Install the Initio® 2 software available at [www.jeulin.fr/Initio@2](http://www.jeulin.fr/Initio@2)
- Connect the Initio® 2 to the computer via the USB- $\mu$ USB cable
- Launch the Initio® 2 software
- Load the saved acquisition by clicking on the icon:

The data then appears on the screen. Right-clicking the chart provides access to various graphic tools.

## 3. Data visualisation

It is possible to view the data at any moment in tabular or chart format. To do so, simply press the Back button (2). The display then toggles between table and chart measurements. Regardless of the current display, data acquisition can be monitored.

After 5 seconds without activity, the table or chart is displayed in full screen. To again obtain the associated functions, turn the knurled knob.

In the table, the arrows situated at the bottom right are used to move between the rows of the table.

## 4. Analogue outputs

Initio® 2 is equipped with banana jack and 3.5 mm jack analogue outputs.

### 4.1 Banana jacks

- Connect the Initio® 2 to an ExAo interface equipped with direct inputs using banana leads.

#### Correspondence table

- **Photometer**

- o Calibre 0 to 2000 lux

<u>-5V</u>	<u>0 lux</u>
<u>5V</u>	<u>2000</u>

- o Calibre 0 to 200klux

<u>-5V</u>	<u>0 lux</u>
<u>5V</u>	<u>200 klux</u>

- **Solarimeter**


<u>-5V</u>	<u>0 <math>W.m^{-2}</math></u>
<u>5V</u>	<u>2000 <math>W.m^{-2}</math></u>

#### Example with the Foxy® console

- Start the generalist

- Place the icon of a direct input on the desired axis
- Select the calibre  $\pm 5V$

- Fill out the customisation tab as below:

- Click on the disk to save the configuration

Configuration is complete, the displayed parameter is now illuminance in lux.

## 4.2 Jack

The 3.5mm jack is to connect the Initio® 2 to a Groomy® interface (ref 181 000) for example. The voltage supplied varies between 0 and 5V.

### Correspondence table

- **Light meter**

- Calibre 0 to 2000 lux

0V	0 lux
5V	1

- Calibre 0 to 200klux

0V	0 lux
5V	200 klux

- **Solarimeter**

0V	0 $W \cdot m^{-2}$
5V	2000 $W \cdot m^{-2}$

## 5. WI-FI communication

Using the WIFI (451 470) module, it is possible to visualise the Initio® 2 data on a tablet or similar.

### 5.1 Installation of the WI-FI module

The WI-FI module is implemented with the Initio® 2 switched off.

- Raise the handle of the Initio® 2
- Remove the back by sliding it down
- Position the WI-FI module at the back
- Push the module firmly into place until it clicks. The handle must be able to retract without forcing.

### 5.2 Viewing on a tablet

Once the WIFI module is installed, the tablet must be configured:

- Press the Menu button
- Select Settings
- Select Network

The module information is displayed on the screen.

- On the tablet, choose the name of the WI-FI network displayed on the Initio® 2.
- Open a web browser
- Enter the address displayed by the Initio® 2 in the address bar of the browser.

## 6. Settings

Various information concerning the Initio® can be obtained by pressing on Menu and selecting Settings → Information

This page displays the serial number of the device and the software version, the battery charge level and the operating language. The latter is changed by turning the knurled knob.

## 7. Update

On our website at the address <http://en.jeulin.fr/en-resources/english-resources/software/free-software.html>, the last update for your Initio® 2 is available. Its version is indicated in the designation.

In order to know the version of your Initio® 2, press the “Menu” button and then “Setup” “Informations”.

## 8. Technical specifications

### Light meter

Calibre	Detail
0 to 2000 lux	±1 %
0 to 200 klux	±1 %

### Solarimeter

Calibre	Detail
0 to 2000 $W.m^{-2}$	±1 %

### Supply:

Li-ion battery with sector adapter and micro-USB cable (supplied)

### Analogue output:

-5 V to +5 V via 4 mm safety banana jacks

0 to +5 V via 3.5 mm jack

### Sizes:

95 x 123 x 44 mm

## 9. After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.

For any repairs, adjustments or spare parts please contact:

### JEULIN – TECHNICAL SUPPORT

468 rue Jacques Monod

CS 21900

27019 EVREUX CEDEX FRANCE

**+33 (0)2 32 29 40 50**