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Test and Measurement

Ref :
281 488


Français – p 1

English – p 4

Version : 8110

Alimentation EVO2[®] R15F5 / 2A
EVO2[®] R15F5 / 2A power supply

1 Description

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
0 825 563 563*

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1 Description

- 1 – Plastic box with double insulation
- 2 – On / off switch
- 3 – Voltage adjustment 0 to +15 V
- 4 – Output sockets 0 to +15 V
- 5 – External control switch
- 6 – Fixed 5 V output

2 Setting up

- Connect the power supply to the 230 V mains supply.
- Set the voltage and current buttons to the left (minimum value)
- Connect your circuit (lamp, motor, voltage divider, etc.) to the selected output sockets.
- Toggle the master switch (2) to " I ". The power-on indicator lamp and the 2 displays light on.
- Adjust voltage and current to desired values

3 Precautions for use

- Avoid using the power supply in direct sunlight or close to a heat source.
- Never block the ventilation vents when the power supply is operating.
- If this device is used in a manner not mentioned in these instructions, the safety ensured by the device may be compromised.

4 Technical characteristic

- **Adjustable output** from 0 to +15 V / 2 A,
- **Fixed output** +5 V / 1 A,
 - Voltage drop less than 20 mV,
 - Residual ripple less than 10 mV,

Note that the GND of the fixed outputs and the GND of the adjustable output (2 black sockets) are electrically separated. They can be linked with a banana lead.

- 2 LCD displays 2.000 points indicate output voltage and current
- Mains voltage : 230 V +/- 10 %, 50-60 Hz,
- Class II double isolation (no ground connexion)
- Consumption : 65 VA maxi,
- Primary protection with 2 A fuse slow blow (not reachable for user),
- Housing dimension: 305 x 305 x 125 mm
- Weight: 4,5 kg

This device complies with:

- The Low Voltage directive 2006/95/CE ;
applied standard : EN 61010-1:2010
- The EMC directive 2004/108/CEE ; applied standards : EN 60326-1

5 Environmental conditions for use of the device

- For indoor use only,
- Altitude less than 2000 meters,
- Operating temperature from 15 °C to 40 °C,
- Maximum relative humidity of 80 % up to 31 °C and 50 % up to 40 °C,
- Mains supply variations of +/- 10 %,
- Installation category II, pollution level 2.

6 Care and maintenance

In case of a malfunction:

- Check the presence of a mains voltage (230 V), the red indicator lamp under the on / off switch is on.
- Check your set-up so as to detect any possible short circuits

If the device still does not work, contact the After Sales Service.

In any cases, **do not try to open the device.**

Store the original packaging; it can be used to bring the device to our workshops for maintenance, if necessary.

7 After-Sales Service

The device is under a 2-year guarantee, it must be sent back to our workshops.

For any repairs, adjustments or spare parts please contact:

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

** Call cost not included.*

*** Only for call from foreign countries.*