ExAO

Capteur Foxy / Tooxy / AirNeXT

Data logging

Sensor

Réf :
482 047


Français – p 1

English – p 4

Version : 8106

Capteur force Foxy
Foxy force sensor

1 Généralités

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
0 825 563 563*

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1 General

A true electronic dynamometer, this sensor allows measuring force, and in particular, investigating weights, moment, concurrent forces, tensile and compression forces.

2 Use

This sensor has to be used with any console in the Foxy range.

2.1 Connection to the Foxy console

- Place the sensor opposite the console connector.
- Slide the connector following the console guide.
- Push home until it clips in.

2.2 Commissioning

The zero must be adjusted with the sensor in position (in the upper or lower position). The positioning depends on the anticipated use.


2.3 Use in stand-alone or connected mode

2.3.1 In stand-alone mode

Adjustment:

Select the calibration then adjust zero:

A long press on the adjustment button: Segment  flashes.

A short press on the adjustment button: Taring OK, segment  goes out.

After adjusting zero, the value displayed is the instantaneous measured value.

2.3.2 When connecting to the computer

The sensor display goes out. The values measured are visible only on the computer screen and the buttons are inactive.

The measuring principle is identical as for stand-alone mode, simply follow the instructions on the computer screen.

3 Additional equipment

No additional equipment is required, except for the attachment systems (and the console).

4 Technical characteristics

- Calibrations: $\pm 5\text{N}$, $\pm 50\text{N}$ (displayed to 4 significant figures)
- Accuracy: 0,1% of full scale
- Resolution : 0,025 N

5 After-Sales Service

The device is under a 2-year guarantee, it must be sent back to our workshops.
For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediately to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.
** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux