Ex.A.O.

Capteur Tooxy / Foxy / AirNeXT

Data logging

Sensor

**Réf :
482 072**

Français – p 1

English – p 4

Version : 8106

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1. Généralité

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
0 825 563 563*

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1. General

This sensor measures the irradiance of UV illumination in two ranges. It also provides the value of the UV index.

2. Use

This sensor functions with Tooxy®, Foxy® and AirNeXT® consoles. There are several possible operating modes:

Standalone mode, the value is displayed directly on the sensor

Mobile mode with the AirNeXT® nomad console

Connected mode when the console is connected to a computer

2.1 Description

This sensor measures irradiance in $W.m^{-2}$ two ranges

- 0 – 150 $W.m^{-2}$
- 0 – 20 $W.m^{-2}$

The change from one range to the other is made by pressing the green button (4)

This sensor also determines the UV index of the received illumination. The change from one function to the other is made by pressing the orange button (2)

2.2 Inserting the sensor

The sensor is connected using the connector on the back of the case:

- align the sensor with the measurement channel (pod or console).
- insert the sensor into the guide system.
- follow the console guide system.
- push in all the way.

2.3 Operation in standalone mode

Start up

Insert the sensor in the measurement channel
Switch on the console or the pod

Display

A value is displayed directly on the sensor. The change from one function to the other or from one sensor to the other is made as described above

Reset

The detection system can be reset.

Close off the light sensitive sensor

Make a long press on the red button (3) until a flashing arrow appears at the bottom left of the display

Make a long press to validate the reset

2.4 Connected to the PC

2.4.1 Foxy® and Tooxy®

Start up

Connect the sensor to the console

Switch the console on

- connect the console to the computer
- once the software is running, choose an application

Parameter settings

Slide the sensor icon on one of the axes

The settings are available in the tabs under the axes.

2.4.2 AirNeXT®

Start up

Insert the sensor in the pod

Connect AirNeXT® to the computer

Switch AirNeXT® and the pod on

Once the software is running, choose an application

Wait until the pod is detected and then validate

Parameter settings

Slide the sensor icon on one of the axes

The settings are available in the tabs under the axes.

Note : In this operating mode, the sensor buttons are deactivated and the settings are made through the software.

2.5 Mobile mode (AirNeXT®)

Start up

Insert the sensor in the pod
Switch AirNeXT® and the pod on
Wait until the pod is detected and then validate

Parameter settings

Select an application
The measured value is displayed at the right of the screen
Press on the value to access sensor parameter settings (functions, range, reset)

3. Other equipment

The sensor is supplied with its light detection system.

4. Technical characteristics

Irradiance function	
Calibres	0 - 20 W/m ² 0 - 150 W/m ²
Resolutions	0,01 W/m ² for 0 - 20 W/m ² 0,1 W/m ² for 0 - 150 W/m ²
UV index function	
Range	1 to 20

5. After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.
For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediately to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.
** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux