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**0 825 563 563\***

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## 1. Principle

The principles of immune reactions are difficult to understand on a purely theoretical level. The handling of models by the students allows a heuristic approach favourable to the acquisition of this complex knowledge.

This model thus allows the student

- *To demonstrate the specificity of antibodies*
- *To simulate the formation of antigen/antibody complexes*
- *To simulate comparatively the vaccinotherapy and serotherapy*

## 2. Composition – set up

The immunology model is composed of pre-cut elements in magnetic material that can be moved and combined on a metalised plate. These elements represent and “materialise” the bacteria and their antigenic patterns on the one hand, the specific antibodies on the other hand.

Using these elements, the student can make several experiments and simulations.

The device contains the following parts:

- 1 support metalised plate A4
- 3 “Bacteria” of type rod-shaped (bacillus)
- 3 “Bacteria” of type Coccus
- 6 Deactivated bacteria antigen
- 3 B lymphocytes type 1
- 3 B lymphocytes type 2
- 6 antibody B type 1
- 6 antibody B type 2
- 1 Storage pouch

## 3. Procedure

Remove the elements from the pre-cut magnetic plate.

Based on the exercises proposed for the students, distribute or ask the students to select the elements necessary for resolving the set problem. The remaining unused parts can be stored in transparent poche.

Some examples of usage are given below.

Based on a simple problem, the student can develop his knowledge regarding the different immunological phenomena.

### 3.1 Identify the elements of the model

The objective is the identification by the students of the different elements represented in the model, with reference to the real that one can study using video images, or of the demonstration model of the teacher, or a scientific text.

These are:

- Antigen and antigen patterns
- Microorganisms (bacteria)
- Antibodies (among others, toxins)

## 3.2 Experiments

### ❖ Objectives

- Understand the different strategies of pathogenic microorganisms
- Show the specificity of antigen patterns
- Show the nature and the specificity of the antigen/antibodies reaction
- Show the principle of help in combating the microbes (vaccinotherapy and serotherapy)

### ❖ Possible activities

- **Simulate an antibody/antigen reaction with agglutination of a variety of bacteria**

This reaction can be studied from an experiment conducted by the student (brucellosis) or presented by the videogram.

After the experiment, the students can prepare diagrams representing the reaction.

- **Simulate an antibody/toxin reaction**

Based on the analysis of a historical text (or representative diagrams), we can study the development of tetanus or diphtheria and the reactions of the organism.

- **Simulate a vaccinotherapy or a serotherapy**

After the analysis of two strategies (diagrams, historical texts, scientific texts...), the student can select from the elements of the model those intervening in each case, then represent the mechanisms involved in a diagram.

## 3.3 Critique of the models / test its limits

Different questions can be posed to the students to help them find the differences between the model and the real and identify the elements of the real that the model does not take into account:

- Relative sizes of elements
- The missing stages
- Production of antibodies
- The fate of agglutinated bacteria
- Recognition of foreign elements
- Other possible types of vaccines (e.g. use of dead microbes)

## 4. After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.  
For any repairs, adjustments or spare parts please contact:

**JEULIN – TECHNICAL SUPPORT**  
468 rue Jacques Monod  
CS 21900  
27019 EVREUX CEDEX FRANCE

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# Direct connection for technical support

A team of experts  
at your disposal  
from Monday to Friday  
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request  
immediatly to provide you  
with the right answers regarding  
your activity field : Biology, Physics,  
Chemistry, Technology.

## Free service\*

**+33 2 32 29 40 50\*\***

\* Call cost not included.

\*\* Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux